SYNGUARD®

Instruction for use

Gloves for single use compliant with PPE regulation (EU) 2016/425, EN 420 and EN ISO 374

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Instructions for use are to be used in combination with the specific product related information on each product packaging. The gloves are sold as a bundled unit within the shipping carton. In case this bundled unit is dismantled and products are sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit.

The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonised European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed after contamination.

EXPLANATION OF THE STANDARDS AND PICTOGRAMS

EN ISO 374-1	Performance levels against permeation are defined on the basis of the breakthrough times (tested according to EN 16523-1:2015) as follows:						
Type A / B / C	Performance level against permeation according to EN ISO 374-1:2016 +A1:2018	1	2	3	4	5	6
	Breakthrough time in minutes	10	30	60	120	240	480
 ABCDEFGHIJKLMNPST	<p>Type A = Breakthrough time of chemical > 30 minutes against at least 6 chemicals on the list. Type B = Breakthrough time of chemical > 30 minutes against at least 3 chemicals on the list. Type C = Breakthrough time of chemical > 10 minutes against at least 1 chemical on the list.</p> <p>TEST CHEMICALS: A = methanol / B = acetone / C = acetonitrile / D = dichloromethane / E = carbon disulfide / F = toluene / G = diethylamine / H = tetrahydrofuran / I = ethyl acetate / J = n-heptane / K = sodium hydroxide 40 % / L = sulfuric acid 96 % / M = nitric acid 65 % / N = acetic acid 99 % / O = ammonia water 25 % / P = hydrogen peroxide 30 % / S = hydrofluoric acid 40 % / T = formaldehyde 37 %</p>						
EN 374-4:2013	The degradation (in %) indicates the change in the puncture resistance of the gloves after exposure to the respective test chemical.						

This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical gloves are suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemicals in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfections.

EN ISO 374-5: 2016  Virus	Tested for permeability according to EN 374-2:2014 Tested for resistance to penetration by blood-borne pathogens according to EN ISO 374-5 / ISO 16604 (virus resistance). Resistance to bacteria and fungi – pass Resistance to virus – pass The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.		Glove if for single-use only and must not be reused.		
			Read the instructions for use carefully before use.		
			Protect from moisture.		Protect from sunlight.
 XXXX	XXXX = Notified Body responsible for EU type examination and supervising ongoing conformity		This is a non-sterile product.		

PRECAUTIONS FOR USE

Before use, always check the gloves for possible mechanical damage, e.g. holes or tears. Do not use damaged gloves. Glove length corresponds to the end use, as the risk to the wrist area is minimal. Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers and other that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. When removing the gloves, avoid allowing the outer surface of the gloves to come in contact with your skin, because the surface may have been contaminated with blood and other bodyfluids. Avoid snapping, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.

INGREDIENTS / HAZARDOUS COMPONENTS

Certain gloves may contain ingredients known to cause skin irritation or allergic reactions in sensitized individuals. Check the warnings on the specific packaging. Formulation will be shared upon request.

TEMPERATURE DELIMITATION

STORAGE INSTRUCTIONS

DISPOSAL INSTRUCTIONS

	The storage area should be cool, dry and dust-free. Protect gloves from UV light sources, sunlight, oxidizing agents and ozone. Store in the original packaging according to the temperature limits indicated on the packaging.	Dispose of according to the valid regulations for the glove material. Gloves contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.
---	---	--

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