Medical Device Regulation Instruction for use


Gloves for single use compliant with MDR regulation (EU) 2017/745 and EN ISO 455

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EN – INSTRUCTION FOR USE Gloves for single use compliant with MDR regulation (EU) 2017/745 and EN ISO 455			
<p>Instructions for use are to be used in combination with the specific product-related information on each product packaging. The gloves are sold as a bundled unit with the shipping carton. In case products in the bundle are unbundled and sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit.</p> <p>The gloves are classified as Medical Device Regulation (MDR) Category I according to MDR Regulation (EU) 2017/745 and have been certified to comply with this regulation based on the applicable harmonized European standards. The gloves meet the EN/ISO standards shown on each specific packaging. This MDR is for single use only and to be disposed after contamination.</p>			
EXPLANATION OF THE STANDARDS AND PICTOGRAMS			
EN 455	European Standard For Single Use Medical Gloves against Blood Penetration and Viruses		
<p>EN 455 Medical Gloves for Single Use covers any glove that could be used for medical work. The glove must adhere to four separate parts before it is considered safe to be used for medical practice.</p> <p>The four parts are:</p> <p>EN 455-1: Requirements and testing of gloves for freedom from holes EN 455-2: Requirements and testing of gloves for physical properties EN 455-3: Requirements and testing of gloves for biological evaluation EN 455-4: Requirements and testing of gloves for shelf life determination</p> <p>Between them they ensure that the glove will be a barrier against micro-organisms, perform effectively without breaking, protect the user from toxic and hazardous materials, and lastly how long a glove will be fit for use.</p>			
	<p>Abbreviation for acceptable quality limit. AQL refers to the quality level that is the worst tolerable process average when a continuing series of lots is submitted for acceptance sampling. At AQL sampling, when the same number of products are sampled, the smaller the value of AQL, the smaller the number of defects allowed.</p>		Gloves for single-use only and non-reusable.
			EU CE Marking for Machinery Directive
			Protect from moisture.
	The CE mark indicates the good's conformity with safety, health, environmental protection, consumer protection, etc. as required by a series of EU directives.		This is a non-sterile product.
PRECAUTIONS FOR USE			INGREDIENTS / HAZARDOUS COMPONENTS
<p>Before use, always check the gloves for possible mechanical damage, e.g. holes or tears. Do not use damaged gloves. Select gloves with a length appropriate for you, as the risk to the wrist area is minimal</p> <p>Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers or others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. When removing the gloves, avoid the contact of gloves' outer surface with your skin, because the surface may have been contaminated with blood or other body fluids. Avoid snagging, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.</p>			<p>Certain gloves may contain ingredients known to cause skin irritation or allergic reactions in sensitive individuals. Check the warnings on the specific packaging. Formulation will be shared upon request.</p>
TEMPERATURE DELIMITATION	STORAGE INSTRUCTIONS	DISPOSAL INSTRUCTIONS	
	<p>The storage area should be cool, dry, and dust-free. Protect gloves from UV light sources, sunlight, oxidizing agents, and ozone. Store in the original packaging according to the temperature limits indicated on the packaging.</p>	<p>Dispose of gloves according to the valid regulations for the glove material. Gloves contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.</p>	

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