SYNGUARD®


Instruction for use

Gloves for single use compliant with PPE regulation (EU) 2016/425, EN ISO 21420 and EN ISO 374

| | | | |
|----|---------------------------|--|-----------|
| EN | Instruction for Use | Gloves for single use compliant with PPE regulation (EU) 2016/425, EN ISO 21420 and EN ISO 374 | Page 3 |
| DE | Gebrauchsanleitung | Einmalhandschuhe konform mit PSA-Verordnung (EU) 2016/425, EN ISO 21420 und EN ISO 374 | Seite 4 |
| FR | Mode d'emploi | Gants jetables conformes à la réglementation sur les EPI (UE) 2016/425, EN ISO 21420 et EN ISO 374 | Page 5 |
| ES | Instrucciones de uso | Guantes de uso único que cumplen con la normativa sobre EPI (UE) 2016/425, EN ISO 21420 y EN ISO 374 | Página 6 |
| IT | Istruzioni per l'uso | Guanti monouso conformi al regolamento DPI (UE) 2016/425, EN ISO 21420 e EN ISO 374 | Pagina 7 |
| RU | Инструкция по применению | Одноразовые перчатки, соответствуют требованиям Регламента ЕС 2016/425 по СИЗ, стандарта EN ISO 21420 и EN ISO 374 | Стр. 8 |
| PT | Instrução de utilização | Luvras para utilização única em conformidade com o regulamento relativo aos equipamentos de proteção individual (UE) 2016/425, EN ISO 21420 e EN ISO 374 | Página 9 |
| DK | Brugsanvisning | Handsker til engangsbrug, der overholder PPE-forordning (EU) 2016/425, EN ISO 21420 og EN ISO 374 | Side 10 |
| NL | Gebruiksaanwijzing | Handschoenen voor eenmalig gebruik conform de PBM-verordening (EU) 2016/425, EN ISO 21420 en EN ISO 374 | Pagina 11 |
| FI | Käyttöohje | Kertakäyttökäsineet, henkilösuojaimista annetun asetuksen (EU) 2016/425, EN ISO 21420:n ja EN ISO 374:n mukaiset | Sivu 12 |
| CZ | Návod k použití | Rukavice pro jednorázové použití v souladu s nařízením o osobních ochranných prostředcích (OOP) (EU) 2016/425, EN ISO 21420 a EN ISO 374. | strana 13 |
| PL | Instrukcja użycia | Rękawice jednorazowe zgodne z rozporządzeniem w sprawie środków ochrony indywidualnej (UE) 2016/425, EN ISO 21420 i EN ISO 374 | Strona 14 |
| SK | Návod na použitie | Rukavice na jedno použitie v súlade s nariadením o osobných ochranných pracovných prostriedkoch (OOP) (EÚ) 2016/425, EN ISO 21420 a EN ISO 374 | strana 15 |
| RO | Instrucțiuni de utilizare | Mănuși de unică folosință conforme cu Regulamentul (UE) 2016/425 privind EIP, EN ISO 21420 și EN ISO 374 | Pagina 16 |

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| EN – INSTRUCTION FOR USE Glove for single use compliant with PPE Regulation (EU) 2016/425, EN ISO 21420 and EN ISO 374 | | | | | | | |
|--|---|--|---|---|--|---|-----|
| <p>Instructions for use are to be used in combination with the specific product-related information on each product packaging. The gloves are sold as a bundled unit with the shipping carton. In case products in the bundle are unbundled and sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit.</p> <p>The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonized European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms, and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed after contamination.</p> | | | | | | | |
| EXPLANATION OF THE STANDARDS AND PICTOGRAMS | | | | | | | |
| EN ISO 374-1 | Performance levels against permeation are defined on the basis of the breakthrough times (tested according to EN 16523-1:2015) as follows: | | | | | | |
| Type A / B / C | Performance level against permeation according to EN ISO 374-1:2016 +A1:2018 | 1 | 2 | 3 | 4 | 5 | 6 |
| | Breakthrough time in minutes | 10 | 30 | 60 | 120 | 240 | 480 |
|  ABCDEFHIJKLMNPST | <p>Type A = Breakthrough time of chemical > 30 minutes against at least 6 chemicals on the list. Type B = Breakthrough time of chemical > 30 minutes against at least 3 chemicals on the list. Type C = Breakthrough time of chemical > 10 minutes against at least 1 chemical on the list. TEST CHEMICALS: A = methanol / B = acetone / C = acetonitrile / D = dichloromethane / E = carbon disulfide / F = toluene / G = diethylamine / H = tetrahydrofuran / I = ethyl acetate / J = n-heptane / K = sodium hydroxide 40 % / L = sulfuric acid 96 % / M = nitric acid 65 % / N = acetic acid 99 % / O = ammonia water 25 % / P = hydrogen peroxide 30 % / S = hydrofluoric acid 40 % / T = formaldehyde 37 %</p> | | | | | | |
| EN 374-4:2013 | The degradation (in %) indicates the change in the puncture resistance of the gloves after exposure to the respective test chemical. | | | | | | |
| <p>This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions using samples taken from the palm only and relates only to the chemical gloves suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemicals in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation resistance can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfection.</p> | | | | | | | |
| EN ISO 374-5: 2016  Virus | <p>Tested for permeability according to EN 374-2:2014 Tested for resistance to penetration by blood-borne pathogens according to EN ISO 374-5 / ASTM F1671 (virus resistance). Resistance to bacteria and fungi – pass Resistance to viruses – pass The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.</p> |  | Gloves for single-use only and non-reusable. | | | | |
| | |  | Read the instructions for use carefully before use. | | | | |
| | |  | Protect from moisture. |  | Protect from sunlight. | | |
|  XXXX | XXXX = Notified Body responsible for EU type examination and supervising ongoing conformity |  | This is a non-sterile product. | | | | |
| Food contact (No more than 30 minutes of contact time and 40 degrees of contact temperature). | | | | | | | |
| PRECAUTIONS FOR USE | | | | | | INGREDIENTS / HAZARDOUS COMPONENTS | |
| <p>Before use, always check the gloves for possible mechanical damage, e.g. holes or tears. Do not use damaged gloves. Select gloves with a length appropriate for you, as the risk to the wrist area is minimal Donning & doffing gloves in a proper way is a skill that needs to be practiced by healthcare workers or others that use gloves. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. When removing the gloves, avoid the contact of gloves' outer surface with your skin, because the surface may have been contaminated with blood or other body fluids. Avoid snagging, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.</p> | | | | | | <p>Certain gloves may contain ingredients known to cause skin irritation or allergic reactions in sensitive individuals. Check the warnings on the specific packaging. Formulation will be shared upon request.</p> | |
| TEMPERATURE DELIMITATION | | STORAGE INSTRUCTIONS | | | DISPOSAL INSTRUCTIONS | | |
|  | | <p>The storage area should be cool, dry, and dust-free. Protect gloves from UV light sources, sunlight, oxidizing agents, and ozone. Store in the original packaging according to the temperature limits indicated on the packaging.</p> | | | <p>Dispose of gloves according to the valid regulations for the glove material. Gloves contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.</p> | | |

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