
# Personal Protective Equipment Instruction for use


Gloves for single use compliant with PPE regulation (EU) 2016/425, EN ISO 21420 and EN ISO 374

EN	Instruction for Use	Gloves for single use compliant with PPE regulation (EU) 2016/425, EN ISO 21420 and EN ISO 374	Page 3
DE	Gebrauchsanleitung	Einmalhandschuhe konform mit PSA-Verordnung (EU) 2016/425, EN ISO 21420 und EN ISO 374	Seite 4
FR	Mode d'emploi	Gants jetables conformes à la réglementation sur les EPI (UE) 2016/425, EN ISO 21420 et EN ISO 374	Page 5
ES	Instrucciones de uso	Guantes de uso único que cumplen con la normativa sobre EPI (UE) 2016/425, EN ISO 21420 y EN ISO 374	Página 6
IT	Istruzioni per l'uso	Guanti monouso conformi al regolamento DPI (UE) 2016/425, EN ISO 21420 e EN ISO 374	Pagina 7
RU	Инструкция по применению	Одноразовые перчатки, соответствуют требованиям Регламента ЕС 2016/425 по СИЗ, стандарта EN ISO 21420 и EN ISO 374	Стр. 8
PT	Instrução de utilização	Luvas para utilização única em conformidade com o regulamento relativo aos equipamentos de proteção individual (UE) 2016/425, EN ISO 21420 e EN ISO 374	Página 9
DK	Brugsanvisning	Handsker til engangsbrug, der overholder PPE-forordning (EU) 2016/425, EN ISO 21420 og EN ISO 374	Side 10
NL	Gebruiksaanwijzing	Handschoenen voor eenmalig gebruik conform de PBM-verordening (EU) 2016/425, EN ISO 21420 en EN ISO 374	Pagina 11
FI	Käyttöohje	Kertakäyttökäsineet, henkilösuojaimista annetun asetuksen (EU) 2016/425, EN ISO 21420:n ja EN ISO 374:n mukaiset	Sivu 12
CZ	Návod k použití	Rukavice pro jednorázové použití v souladu s nařízením o osobních ochranných prostředcích (OOP) (EU) 2016/425, EN ISO 21420 a EN ISO 374.	strana 13
PL	Instrukcja użycia	Rękawice jednorazowe zgodne z rozporządzeniem w sprawie środków ochrony indywidualnej (UE) 2016/425, EN ISO 21420 i EN ISO 374	Strona 14
SK	Návod na použitie	Rukavice na jedno použitie v súlade s nariadením o osobných ochranných pracovných prostriedkoch (OOP) (EÚ) 2016/425, EN ISO 21420 a EN ISO 374	strana 15
RO	Instrucțiuni de utilizare	Mănuși de unică folosință conforme cu Regulamentul (UE) 2016/425 privind EIP, EN ISO 21420 și EN ISO 374	Pagina 16

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EN – INSTRUCTION FOR USE Glove for single use compliant with PPE Regulation (EU) 2016/425, EN ISO 21420 and EN ISO 374							
<p>Instructions for use are to be used in combination with the specific product-related information on each product packaging. The gloves are sold as a bundled unit with the shipping carton. In case products in the bundle are unbundled and sold separately, the distributor must ensure that the instructions for use are accompanied with each separate unit.</p> <p>The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonized European standards. These gloves are designed to provide protection against specific chemicals tested, micro-organisms, and particulate radioactive contamination (if applicable). The gloves meet the EN/ISO standards shown on each specific packaging. This PPE is single-use only and to be disposed after contamination.</p>							
EXPLANATION OF THE STANDARDS AND PICTOGRAMS							
EN ISO 374-1	Performance levels against permeation are defined on the basis of the breakthrough times (tested according to EN 16523-1:2015) as follows:						
Type A / B / C	Performance level against permeation according to EN ISO 374-1:2016 +A1:2018	1	2	3	4	5	6
	Breakthrough time in minutes	10	30	60	120	240	480
 ABCDEFGHIJKLMNOST	<p>Type A = Breakthrough time of chemical &gt; 30 minutes against at least 6 chemicals on the list.            Type B = Breakthrough time of chemical &gt; 30 minutes against at least 3 chemicals on the list.            Type C = Breakthrough time of chemical &gt; 10 minutes against at least 1 chemical on the list.  <b>TEST CHEMICALS:</b>            A = methanol / B = acetone / C = acetonitrile / D = dichloromethane / E = carbon disulfide / F = toluene / G = diethylamine / H = tetrahydrofuran / I = ethyl acetate / J = n-heptane / K = sodium hydroxide 40 % / L = sulfuric acid 96 % / M = nitric acid 65 % / N = acetic acid 99 % / O = ammonia water 25 % / P = hydrogen peroxide 30 % / S = hydrofluoric acid 40 % / T = formaldehyde 37 %</p>						
EN 374-4:2013	The degradation (in %) indicates the change in the puncture resistance of the gloves after exposure to the respective test chemical.						
<p>This information does not reflect the actual duration of protection at the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance have been assessed under laboratory conditions using samples taken from the palm only and relates only to the chemical gloves suitable for the intended use because the conditions (such as temperature, abrasion, and degradation) at the workplace may differ from the testing conditions. Used gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemicals in contact with the gloves, etc. may shorten the actual service life of gloves significantly. For corrosive chemicals, degradation resistance can be the most important factor to consider when selecting chemical-resistant gloves. Before use, inspect the gloves for any defect or imperfection.</p>							
EN ISO 374-5: 2016  Virus	<p>Tested for permeability according to EN 374-2:2014            Tested for resistance to penetration by blood-borne pathogens according to EN ISO 374-5 / ASTM F1671 (virus resistance).            Resistance to bacteria and fungi – pass            Resistance to viruses – pass            The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.</p>		Gloves for single-use only and non-reusable.				
			Read the instructions for use carefully before use.				
			Protect from moisture.		Protect from sunlight.		
 XXXX	XXXX = Notified Body responsible for EU type examination and supervising ongoing conformity		This is a non-sterile product.				
PRECAUTIONS FOR USE						INGREDIENTS / HAZARDOUS COMPONENTS	
<p>Before use, always check the gloves for possible mechanical damage, e.g. holes or tears. Do not use damaged gloves. Select gloves with a length appropriate for you, as the risk to the wrist area is minimal. Donning must be performed in the correct order to prevent transmission of infections. Keep hands clean before donning gloves. When removing gloves, avoid the contact of gloves' outer surface with your skin, because the surface may have been contaminated with blood or other body fluids. Avoid snagging, as this may cause contaminants to splash into your eyes or mouth or onto your skin or other people nearby.</p>						<p>Certain gloves may contain ingredients known to cause skin irritation or allergic reactions in sensitive individuals. Check the warnings on the specific packaging. Formulation will be shared upon request.</p>	
TEMPERATURE DELIMITATION	STORAGE INSTRUCTIONS		DISPOSAL INSTRUCTIONS				
	<p>The storage area should be cool, dry, and dust-free. Protect gloves from UV light sources, sunlight, oxidizing agents, and ozone. Store in the original packaging according to the temperature limits indicated on the packaging.</p>		<p>Dispose of gloves according to the valid regulations for the glove material. Gloves contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.</p>				

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