

General Isolation Gown - 45GSM

Features

- Made from lightweight polypropylene and multi-layer materials
- Neck tie and elastic wrists
- Suitable for isolation and basic protection from bacteria and particulate
- Strong sealed seams to withstand tough conditions
- Made of a clear laminate material
- Breathable material for extended wear and fluid-proof barrier protection



Applications

- Medical Purpose / Examination
- Industrial purpose / PPE
- Laboratory
- Healthcare and nursing
- General housekeeping
- IT Industry



SKU	Size	Color	Material	Package	Carton Dimension
IG45SN1024	115*137cm	White	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG45SN1025	120*140cm	White	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG45SN1014	115*137cm	Blue	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG45SN1015	120*140cm	Blue	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG45SN1064	115*137cm	Yellow	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG45SN1065	120*140cm	Yellow	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm



Round Collar



Collar Tie



Belt



Elastic Cuff