

Isolation Gown - AAMI Level 2

Features

- Made from lightweight SMMS fabrics
- Neck tie and knitted cuffs
- Suitable for isolation and basic protection from bacteria and particulate AAMI level 2 performance
- Ultrasonic welding between sleeves and body
- Blue and Yellow color available

Quality standards

- ANSI/AAMI PB70-2012 :Level 2
- AATCC 42-2017



Applications

Medical Purpose / Examination

Industrial purpose / PPE

Laboratory

Healthcare and nursing

General housekeeping

IT Industry



SKU	Size	Color	Material	Package	Carton Dimension
IG25U21024	115*137cm	White	SMMS 25GSM	100pcs/bag,10bags/ctn	610*260*300mm
IG25U21025	120*140cm	White	SMMS 25GSM	100pcs/bag,10bags/ctn	610*260*300mm
IG25U21014	115*137cm	Blue	SMMS 25GSM	100pcs/bag,10bags/ctn	610*260*300mm
IG25U21015	120*140cm	Blue	SMMS 25GSM	100pcs/bag,10bags/ctn	610*260*300mm
IG25U21064	115*137cm	Yellow	SMMS 25GSM	100pcs/bag,10bags/ctn	610*260*300mm
IG25U21065	120*140cm	Yellow	SMMS 25GSM	100pcs/bag,10bags/ctn	610*260*300mm



Armhole



Collar Tie



Belt



Knitted Cuff