


Isolation Gown - AAMI Level 3

Features

- Comfortable and waterproof materials with classic tie waist
- Stitching seam sealed with waterproof tapes
- Suitable for isolation and basic protection from bacteria and particulate AAMI Level 3 performance
- Neck tie and elastic wrists
- Blue, White and Yellow color available



Quality Standards

- ANSI/AAMI PB70:2012
- AATCC 42-2017
- AATCC 127-2018

Primary Industry

- Moderate risk to be used, for example, during arterial blood draw, inserting an Intravenous (IV) line, in the Emergency Room, or for trauma cases
- Non-Surgical use

Applications

- Medical Purpose / Examination
- Industrial purpose / PPE
- Laboratory
- Healthcare and nursing
- General housekeeping
- IT Industry



SKU	Size	Color	Material	Package	Carton Dimension
IG4531015	120*140cm	Blue	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG4531025	120*140cm	White	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm
IG4531065	120*140cm	Yellow	SPP+PE 45GSM	100pcs/bag,10bags/ctn	500*300*450mm



Armhole



Belt



Knitted Cuff