

06 DISPOSABLE HEAT PATCH

Using heat patch is the most convenient and safest way to provide skin-touching continuous and stable warm. It is portable and has great effect in preventing from frostbite in winter. Just one single piece can resist full day's cold.





HEAT PATCH



Main features

It's disposable and easy to provide relief of minor muscular and joint aches and pains associated with over exertion, strains, sprains and arthritis, period pain, etc.

It will gradually get warm and last for at least 6/8 hours.

Instruction

Open the outer pouch, leaving the internal pack to react with air. It needs some time to reach its therapeutic temperature. Please refer to the instruction for use for specific product usage.

				
Item Number	L095130A16	L055950A8	L070095A8	H080185B0.5
Size	Cover:105.5*15.5CM Heat patch:13*9.5CM	9.5*5.5CM	9.5*7CM	8*18.5CM
Weight	38±5G	30G	21G	11±2G
Material	Non-woven, PP Non-woven, Heating material	PP Non-woven, Heating material	PP Non-woven, Heating material, PE stick on clothes/PP film	Non-woven, Heating material
Using Position	Waist	Hand	Pelma	Eye

				
Item Number	H055095B8/L095130A12	L100150A8	H095285B8	H095285B8/H086260B8/H085250B8
Size	9.5*5.5CM/9.5*13CM	15*10CM	28.5*9.5CM	28.5*9.5/26.5*8.6/25*8.5CM
Weight	30±5G/38±5G	30G	55±5G	55±5G/51±5G/46±5G
Material	Spunlace PE/PP Non-woven, Heating material, PE Stick on skin/Spunlace stick on skin/ PE stick on clothes/PP film			PP Non-woven, Heating material, Ivory board film, PE stick on clothes
Using Position	Multi-use	Belly	Shoulder/Back	Foot